

FLOW CHART

Instant Essential Amino Acid(EAA)-M



L-LEUCINE,L-ISOLEUCINE,L-VALINE,L-ALANINE,L-ARGININE,L-ASPARTIC ACID,L-CYSTINE,L-GLUTAMIC ACID,L-GLYCINE,L-HISTIDINE,L- LYSINE HCL,L-METHIONINE,L-PHENYLALANINE,L-PROLINE,L-SERINE,L-THREONINE,L-TRYPTOPHAN,L-TYROSINE

Mixing (Normal Temperature)

Microwave Drying(80°C)

Sieving

Mixing

Packing

Storage